BSPS COSHH Risk assessment

POCT-BSPS-COSHH-001		
Site: RBH, ASPH, RSH, WPH and FPH	Department: POCT	
Name of substance (this should include the trade name and chemical name):	Physical form (aerosol; gas; fumes; dust; liquid; powder etc)	
Radiometer SDS M1 038 sensor cassettes:	Cassette containing inner liquid electrolyte and salts.	
946-005 SC90 300 BG, LYT, MET, OXI + QC		
2 946-009 SC90 900 BG, LYT, MET, OXI + QC		
946-060 SC90 1200 BG,LYT,MET,OXI+QC		
946-010 SC90 100 BG, LYT, MET, OXI +QC		
946-008 SC90 600 BG, LYT, MET, OXI + QC		
946-705 SC90 Ki 300 BG, LYT, MET, UREA, CREA, OXI + QC		
CAS Registration number:	Additional comments:	
None		
Suppliers name:	Suppliers contact information:	
Radiometer Ltd	Manor Court	
	Manor Royal, Crawley	
	Crawly, West Sussex RH10 9FY	
	Tel: +44 (0) 1293 517599	
	www.radiometer.co.uk	
Description of the process where chemical/substance/product is used:		
Sensor cassettes are self-contained ABL90 FP packs for the measurement of the compartment inside the machine. The machines are located in acute, medical, s	e parameters shown on the label of the box that contains the sensor cassette. They are placed in a surgical and urgent care wards across all BSPS POCT sites.	

Is there a safe system of work for the activity?	Is it possible for the hazardous substance to be substituted for a safer alternative? (Yes (No)
--	---

Radiometer ABL90 Flex Plus Sensor Cassettes

POCT-BSPS-COSHH-001 | Version: 1.0 | Authorised By: David Baulch | On 17-Nov-2022

LCOSHH Risk Assessment form GEN-BSPS-H&S-2 | Version: 3.1 | Authorised By: Elaine Inglis | On 27-Sep-2022

Document ID: 401-118901125 | Printed: 10-Nov-2023 at 16:41

Safety/Risk Phrases	Hazard Classification (corrosive; environmental hazard; irritant; oxidising agent etc)	Route of entry (absorption; ingestion; inhalation; injection etc)	On each occasion quantity used:	On each occasion duration of exposure (minutes)	Frequency of exposure
None	Under the REACH regulation, this product is defined as an article with no intended release of substances. The product is not covered by the regulations on classification and labelling of chemical substances and materials.	Absorption and ingestion but not relevant due to design of product unless cassette breaks. Skin and eye contact by absorption. Ingestion not likely, due to the form of the product.	Small quantities. Cassette size is 2ml	Staff will not be in direct contact with contents as they are self-contained unless damage to cassette occurs. Sample processing takes less than a minute and maintenance varies but is less than 30minutes. During sample processing cassette is locked in its compartment at all time.	Daily. Due to the design of the product, exposure is likely to be brief and in small quantities.
Individuals/Groups exposed ?	Patients/Staff/Visitors/Contractor	Numbers of individuals potentially of	exposed. 40 POCT staff	who do maintenance and en	gineers.
Does the substance have a Workplace Exposure Limit (WEL)? Yes No If there is a WEL, what is it? (EH40 reference information to be included i.e.)					
irritate Ingestion: Not likely du	ct, exposure is likely to be brief and in small quan	utities. Skin contact: Prolonged or repe	ated contact may cause	irritation. Eye contact: Direc	t contact may
		Controls			
Describe the existing controls i	in place (your description should include ventilati estrict access; staff rotation etc)	ion such natural, mechanical and LEVs;	segregation; substitutio	on to less hazardous substand	ce; use of
	ct, exposure is likely to be brief and in small quan oid contact with skin and eyes. Cassette changes aled to avoid spillages if any.				l in SULO bins in

is required and the current arrangements with Trust Occupational Health	Atmospheric monitoring: Is atmospheric monitoring required? If yes, please outlined the process currently in use. Please ensure that all records of atmospheric monitoring are retained in line with recommended record retention period.		
Not required	Not required		
Storage and transport:	Disposal:		
and store in POCT fridge until required.	Disposed in UN3291 clinical waste SULO bins with a lid on in its original container. Any broken cassettes are collected in containers and sealed securely. Due to the design of the product, exposure is likely to be brief and in small quantities		
PPE used:			
	Image: Dust MaskImage: Dust MaskImage		
	Emergency procedures		
First Aid:	Management of Spillages		
Inhalation: Not relevant	Absorb spillage with suitable absorbing material. Flush away small spillages with plenty of water. Personal precautions: Avoid contact with skin and eyes. Environmental precautions: Avoid discharge into drains, water courses or onto		

Legal entity host: Frimley Health NHS Foundation Trust

Initial Risk Level Assessment Is the control of the risk adequate? It is important to include more than one level of risk if the substance is used in multiple processes. Please use the risk matrix appropriate for the site you are located at.

Consequence/impact X Likelihood/probability = **Insert** numerical value and grade (low, medium, high) Unlikely. Risk is 2-low risk

Where risk level is identified as being high following implementation of controls has this been escalated to senior management and discussed in BSPS Health and safety meetings Yes/No

Action plan (what additional actions are to be taken if the risk level is identified as high or very high)

Proposed action required	By Whom please include title of manager who has been informed	Action Start Date	Action Due Date	Date completed

Risk Level Assessment following additional actions:

Is the control of the risk now deemed as adequate? It is important to include more than one level of risk if the substance is used in multiple processes. Please use the risk matrix appropriate for the site you are located at.

Consequence/impact X Likelihood/probability = Insert numerical value and grade (low, medium, high)

Are the risks now adequately controlled? Yes/No

Communication and staff training:

Induction, training, competency, SOPs, iPassport skill.

Assessor Name:	Assessor Title:	Date:	Review date:
Eulitty Gumunyu	POCT Coordinator and H&S Officer	07/11/22	21/11/23
Managers Name:	Managers Title	Date:	Review date:
David Baulch	POCT Manager	16/11/23	21/11/23
	_		

Radiometer ABL90 Flex Plus Sensor Cassettes

POCT-BSPS-COSHH-001 | Version: 1.0 | Authorised By: David Baulch | On 17-Nov-2022

Legal entity host: Frimley Health NHS Foundation Trust